



Mortuary Fridge Band



INFORMATION SHEET

For use by mortuary technicians to secure arms and legs of deceased within the confines of their fridge tray.



Mortuary Fridge Band

The Hospital Direct Mortuary Fridge Band was specifically designed so that mortuary technicians can secure arms and legs of deceased within the confines of their fridge tray.

WHY THEY WERE DEVELOPED

Fitting the arms and legs of large and bariatric deceased onto a fridge tray has long caused problems for mortuary technicians as they can expand over the sides of the trays. To this end staff have had to develop their own, not always, dignified solutions and there was a definite need for a product that would solve the problem quickly and easily whilst also ensuring the dignity of the deceased.

OUR SOLUTION

The Hospital Direct Mortuary Fridge Band is a unique multi adjustable product enabling mortuary staff to carefully secure the arms or legs of a medium to large deceased patient within the confines of their fridge tray. They can then be safely transferred in and out of their fridge compartment knowing that the risk of body contact around the fridge during this manoeuvre has been dramatically reduced.

Quick and easy to use it slides directly under a body bag or wrapped body and reduces risk of poor handling and over reaching. Multiple bands can be used.


PRODUCT DETAILS

This specialist band is 200mm wide and adjustable in length so that it can fit all deceased. It has a 2 point adjustable buckle system allowing adjustment from small to bariatric. For example: it can be closed up to 510mm to go around smaller patients and opened out to 760mm for larger patients.

British made using specially chosen laminated materials, strong nylon webbing and large plastic buckles for comfortable adjustment. **Buckles tested to -40°C.**

It can be laundered up to 74°C and tumble dried at a low temperature.

Disinfectant wipes can also be used in accordance with your organisation's cleaning protocol and chlorine bleach is not to be used.



| Mortuary Fridge Band | Product Code | Size | Overall Size |
|----------------------|--------------|-----------------------|------------------------|
| Mortuary Fridge Band | MFB001 | One size (adjustable) | 200mm (w) x 1450mm (l) |

INSTRUCTIONS FOR USE

1

Unfasten the black buckle if it's connected so the band lays flat.

2

Slide the band under the body from the head or foot end down to the area that needs supporting.

3

Before connecting the band, make sure there's an equal length both sides.

4

Bring together both halves of the band round the body and fasten securely using one of the two buckles which best fits that particular body.

5

If the band is still loose further adjustment can be made by pulling the black webbing strap through the buckle until the band is securely fastened.

6

Once the band is finished with unfasten the buckle and slide the product out from under the body.

Visit our Website
www.hospitaldirect.co.uk

Wash at 74°C
Tumble dry on low heat
Iron on low heat
Do not use bleach
Do Not Use fabric conditioner



+44(0)1588 640144

3

RELATED PRODUCTS

XCcube® V3

New for 2021 the XCcube® V3 ensures that your deceased patient is always treated with respect and dignity. This unique expandable cube offers the best end of life care in terms of dignity and discretion and revolutionises the way staff move deceased to mortuaries. This newest version improves on the last with new features, greater ease of use, an improved transport trolley and increased robustness.




Mortuary Webbing Slings

The Mortuary Webbing Slings have been specifically designed so that mortuary technicians can safely lift and lower a bariatric patient into a coffin using a five point overhead flat lifting gantry system. Each sling is made from a heavy duty 90mm wide simplex polyester webbing and has a SWL of 1.5 tonnes. They are CE marked and conform to BSEN 1492-1:2000.



Surgical Gowns

UK made Disposable Non-sterile Surgical Gowns. Made in our premises in Shropshire, this superb range of non-sterile isolation gowns are alcohol repellent, have set in sleeves, adjustable neck sizes and welded seams in areas most at risk of exposure. Available in boxes of 25 in blue or green. We offer four generous sizes from Small to XLarge. CE marked ISO EN 13795.



Hospital Direct is a Shropshire based UK manufacturer and supplier of a wide range of Patient Handling Equipment, PPE and other innovative solutions for healthcare. We have been supplying the NHS and healthcare customers for over 25 years and have an enviable reputation for quality, delivery, customer care and service. We supply across the UK and export around the world. If you have any questions or need advice please either call or email where we will be delighted to help.

Packaging: To reduce wastage we use recyclable cardboard boxes for shipping. If you require different packaging, please call us to discuss.

Disposal: Dispose of in the normal way in accordance with your organisational procedures. If they have been used on highly infected patients please dispose of as such.



Contact our sales team

Our customer service and sales teams are always available to answer any questions you may have. Please call or email us and we will endeavour to answer as soon as possible.

Visit our website / Get a quote

www.hospitaldirect.co.uk

Email

For sales orders: sales@hdigroup.co.uk

For information: info@hdigroup.co.uk

Phone

+44 (0)1588 640144

Address

Hospital Direct (Marketing) Ltd
Units 3 & 4, The Green
Clun, Shropshire SY7 8LG



BAB logos are applicable to the company management systems only and not the product