



EXPERTS IN PATIENT HANDLING SOLUTIONS



PATIENT TRANSFER SHEET – WASHABLE

This guide covers all products with the prefix TSS in the product code.


The **Patient Transfer Sheet – Washable** enables users to transfer their patients from one lying surface to another across a lateral transfer board. Available with 6 or 8 handles.

Strong silicone coated nylon fabric



Polypropylene handles
(available with 6 or 8 handles)

USER INSTRUCTIONS



1

The transfer must be well planned and all staff briefed. Check brakes, height and position before commencement.
Position the **Patient Transfer Sheet** under the patient. Align both surfaces prior to transfer and secure brakes.

2

Using the **Patient Transfer Sheet** to help, log roll the patient and insert a lateral transfer board sufficiently under the patient between the surfaces, to ensure a smooth transfer, making sure that any gap between the surfaces is covered.

3

Depending on the number of staff, select the most appropriate handles to ensure an effortless transfer. Only transfer the patient on safe, co-ordinated commands.

4

When Transfer is complete, remove the sheet by log rolling the patient.

MATERIAL CONTENT

Fabric: 100% Silicone coated nylon

Handles: Polypropylene



STORAGE AND HANDLING

- Store in a clean dry environment that is not subject to excessive heat.
- Keep away from sunlight.
- Store flat and avoid placing heavy objects on top of the product.

WARNINGS AND PRECAUTIONS

- Do not leave the on the floor when not being used.
- Inspect the product, including handles, before use for signs of damage or tears.
- Safe Working Load (SWL): 127kgs



Keep dry



Keep away from sunlight



Caution



Do not use if package is damaged



Storage Temperature Limit

CLEANING INSTRUCTIONS

- Wash as per laundry instructions.
- Can be wiped clean, between use with the same patient, with natural water and detergent or wipes that don't contain alcohol, solvents, bleaching or abrasive agents.
- Cleaning materials used should be patient safe and biodegradable..



Wash at 30-74°C.



Tumble dry on low heat.



Do not use fabric conditioner

DISPOSAL INSTRUCTIONS

- Non-recyclable.
- Dispose of as clinical waste if it has been used by an infectious person.
- Otherwise dispose of through normal waste management.



MANUFACTURER

Hospital Direct (Marketing) Ltd, Units 2, 3 & 4, The Green, Clun, Shropshire SY7 8LG, United Kingdom

EU REPRESENTATIVE: Obelis S.A., Bd Général Wahis 53, 1030 Brussels, Belgium



www.hospitaldirect.co.uk


HOSPITAL DIRECT (HD) ADDITIONAL SYMBOLS

At HD we understand that not everyone is familiar with all equipment, so to help therapists and healthcare professionals in assessing a product's suitability for a client, we have included two assessment tools, Functional Independence Measure (FIM) and Easy Guide Symbols. These will enable the healthcare professional to assess their client and decide whether or not a product will help the person based on their level of dependence and functional ability.

FUNCTIONAL INDEPENDENCE MEASURE (FIM) ASSESSMENT

Our FIM assessment guide for each product will help you decide the suitability of this product for the person's ability and need. Based on the standard criteria from Level 6 where the person can use the product unaided and unsupervised to Level 1 where all the assistance is provided by the carer and the client can do nothing, this guide is easy, quick and simple to use to check suitability against ability and circumstance.

FIM LEVELS QUICK GUIDE. SUBJECT VS CARER HELP



©2018 Hospital Direct (Marketing) Ltd

EASY GUIDE SYMBOLS

We use four symbols to indicate where a product is suitable and safe for the person to use unassisted, with a single carer or with multiple carers. These symbols indicate the minimum recommendation.



HD has a policy of continuous development and, as such, reserves the right to alter specifications (including measurements, materials and colours) without prior notice. If any serious incident has occurred in relation to the device it should be reported to the manufacturer and the competent authority in which the user and/or patient is established.