

EXPERTS IN PATIENT HANDLING SOLUTIONS



ESSENTIAL GLIDE AND LOCK RANGE


This guide covers GL4045E GL4545E GL4650E, GL6580E, GL4045EW, GL4650EW, GL4543SF and GL6580SF.

An economy range of comfortable, hard wearing products that enable the user to glide to the back of the chair and then lock securely in position, preventing sliding forward.

PLEASE NOTE: The instructions below features our **Glide and Lock - Standard** but the **Essential Glide and Lock Range** is used in the same way.

USER INSTRUCTIONS





1 UNASSISTED

Open up your **Essential Glide and Lock** and hold it with your hand through the loop, ensuring the anti-slip panel is on the top of the product **ON THE REVERSE**.
The directional label will be on the left hand side of the product as you look at it.

Place the **Essential Glide and Lock** on the chair or wheelchair prior to seating the patient/client. Check that the direction of the inner locking strips is correct (see directional arrows on the label).

The normal positioning is for the anti-slip panel to be face down against the chair surface, whilst the label is on the top surface and on the carer's left hand side.





2


Ask the client to sit on the Glide and Lock, making sure that their bottom does not overshoot the product. The best position is for the user to sit forward with their nose over their toes with their hands on the chair or wheelchair arms.
For those people who need a little more grip with their feet, our **Anti-slip Foot Pad** or **HD Deluxe™ Footstool** are extremely useful.



By using their hands and feet to push back, the user will be able to glide to the back of the chair. The locking strips inside the **Essential Glide and Lock** will now lock them into place and prevent them from slipping forwards.

This unassisted method is only suitable for those with sufficient strength in their lower limbs and arms.





3 WITH ASSISTANCE

If the user needs a little more help, you can try the method below, having carried out a risk assessment to make sure that this method is suitable for the patient or client.
When the user is sitting with 'nose over toes', the carer can use a pillow placed just below the knees to push and help glide them back.

The user will now be at the back of the chair or wheelchair and the inner locking strips will lock them into place preventing them sliding forward.

For those who need more assistance, try our **Essential Glide and Lock with Handles**.



VARIETIES AVAILABLE WITHIN THE ESSENTIAL GLIDE AND LOCK RANGE

Breathable

The breathable, quilted, polycotton top fabric layer makes it robust, hard-wearing and comfortable.



Spacer Fabric

The spacer fabric used in this range offers a little extra pressure care and hence comfort for its user.



Wipeclean

The wipeclean top layer fabric makes this more suitable for those users who might be incontinent.



MATERIAL CONTENT

Fabric: **Breathable:** Polycotton/Polyester **Spacer Fabric:** 100% Polyester **Wipeclean:** PVC coated Polyester
Inner Locking Strips: 100% Nylon
Anti-slip panel: PVC coated polyester, phthalate free

STORAGE AND HANDLING

- Store in a clean dry environment that is not subject to excessive heat.
- Keep away from sunlight.
- Store flat and avoid placing heavy objects on top of the product.

WARNINGS AND PRECAUTIONS

- Do not leave or drop on the floor; it is a slip hazard.
- Test for slippiness before use. If over-laundered and doesn't slide please discard or replace.
- Inspect product before use for signs of damage or tears.

CLEANING INSTRUCTIONS

- Wash as per laundry instructions.
- Can be wiped clean, between use with the same patient, with natural water and detergent or wipes that don't contain alcohol, solvents, bleaching or abrasive agents.
- Cleaning materials used should be patient safe and biodegradable.

DISPOSAL INSTRUCTIONS

- Non-recyclable.
- Dispose of as clinical waste if it has been used by an infectious person.
- Otherwise dispose of through normal waste management.



Keep dry



Keep away from sunlight



Caution



Do not use if package is damaged



Storage Temperature Limit



Wash at 30-74°C. Tumble dry on low heat. Do not Iron.
Do not use bleach. Do not use fabric conditioner.



Patient Information website



Medical Device



MANUFACTURER

Hospital Direct (Marketing) Ltd, Units 2, 3 & 4, The Green, Clun, Shropshire SY7 8LG, United Kingdom

EU REPRESENTATIVE: Obelis S.A., Bd Général Wahis 53, 1030 Brussels, Belgium



www.hospitaldirect.co.uk

HOSPITAL DIRECT (HD) ADDITIONAL SYMBOLS

At HD we understand that not everyone is familiar with all equipment, so to help therapists and healthcare professionals in assessing a product's suitability for a client, we have included two assessment tools, Functional Independence Measure (FIM) and Easy Guide Symbols. These will enable the healthcare professional to assess their client and decide whether or not a product will help the patient based on their level of dependence and functional ability.

FUNCTIONAL INDEPENDENCE MEASURE (FIM) ASSESSMENT

Our FIM assessment guide for each product will help you decide the suitability of this product for the person's ability and need. Based on the standard criteria from Level 6 where the person can use the product unaided and unsupervised to Level 1 where all the assistance is provided by the carer and the client can do nothing, this guide is easy, quick and simple to use to check suitability against ability and circumstance.

FIM LEVELS QUICK GUIDE. SUBJECT VS CARER HELP



©2018 Hospital Direct (Marketing) Ltd

EASY GUIDE SYMBOLS

We use four symbols to indicate where a product is suitable and safe for the person to use unassisted, with a single carer or with multiple carers. These symbols indicate the minimum recommendation.



Person Unassisted



Single Carer



Two Carers



Four Carers

HD has a policy of continuous development and, as such, reserves the right to alter specifications (including measurements, materials and colours) without prior notice. If any serious incident has occurred in relation to the device it should be reported to the manufacturer and the competent authority in which the user and/or patient is established.