



EXPERTS IN PATIENT HANDLING SOLUTIONS



HDDELUXE™

HD DELUXE™ – COMFORT SUPPORT PILLOW

This guide covers the product code CLT.

Our **HD Deluxe™ – Comfort Support Pillow** Is designed to provide your end users with the ultimate in head and neck support, comfort and softness. The memory foam used is deliberately chipped so as to allow more flexibility than traditional memory foam pillows. In these pillows, each particle of memory foam reacts directly to the users head shape and weight , creating a super soft surface to sleep on that is specific to them and that can mould to support their particular needs, providing pressure relief at all times. The covers are fully welded, waterproof, vapour permeable. They also reduces the risk of pressure sores developing as the particles are always moving.

USER INSTRUCTIONS



Covered with a soft 4-way stretch fabric designed specifically for direct skin contact, our Comfort Support Pillow is filled with chipped viscoelastic memory foam and conforms to the users head shape for ultimate comfort and support and helps prevent pressure sores developing.



All seams are fully welded for infection control.



Waterproof, the covers are also vapour permeable, biocompatible and antibacterial.

MATERIAL CONTENT

Cover Fabric: 100% polyester
Cover Coating: Polyurethane
Chipped Viscoelastic Foam: Polyurethane
Biocompatible cover (Oekotex certified) can be used directly on the skin



STORAGE AND HANDLING

- Store in a clean dry environment that is not subject to excessive heat.
- Keep away from sunlight.
- Store flat and avoid placing heavy objects on top of the product.

WARNINGS AND PRECAUTIONS

- Do not leave the on the floor when not being used.
- Do not use on a wet floor.
- Always test product grip on your chosen surface before using.
- Inspect before use for signs of wear or damage.



Keep dry




Keep away from sunlight



Caution



Do not use if package is damaged



Storage Temperature Limit

CLEANING INSTRUCTIONS

- Wipeclean with natural water and detergent or wipes that don't contain alcohol, solvents, bleaching or abrasive agents.
- Cleaning materials used should be patient safe and biodegradable.



Wipeclean Only

DISPOSAL INSTRUCTIONS

- Non-recyclable.
- Dispose of as clinical waste if it has been used by an infectious person.
- Otherwise dispose of through normal waste management.



Medical Device



MANUFACTURER

Hospital Direct (Marketing) Ltd, Units 2, 3 & 4, The Green, Clun, Shropshire SY7 8LG, United Kingdom



EU REPRESENTATIVE: Obelis S.A., Bd Général Wahis 53, 1030 Brussels, Belgium



www.hospitaldirect.co.uk

HOSPITAL DIRECT (HD) ADDITIONAL SYMBOLS

At HD we understand that not everyone is familiar with all equipment, so to help therapists and healthcare professionals in assessing a product's suitability for a client, we have included two assessment tools, Functional Independence Measure (FIM) and Easy Guide Symbols. These will enable the healthcare professional to assess their client and decide whether or not a product will help the person based on their level of dependence and functional ability.

FUNCTIONAL INDEPENDENCE MEASURE (FIM) ASSESSMENT

Our FIM assessment guide for each product will help you decide the suitability of this product for the person's ability and need. Based on the standard criteria from Level 6 where the person can use the product unaided and unsupervised to Level 1 where all the assistance is provided by the carer and the client can do nothing, this guide is easy, quick and simple to use to check suitability against ability and circumstance.

FIM LEVELS QUICK GUIDE. SUBJECT VS CARER HELP



EASY GUIDE SYMBOLS

We use four symbols to indicate where a product is suitable and safe for the person to use unassisted, with a single carer or with multiple carers. These symbols indicate the minimum recommendation.



Person Unassisted



Single Carer



Two Carers



Four Carers

HD has a policy of continuous development and, as such, reserves the right to alter specifications (including measurements, materials and colours) without prior notice. If any serious incident has occurred in relation to the device it should be reported to the manufacturer and the competent authority in which the user and/or patient is established.