



EXPERTS IN PATIENT HANDLING SOLUTIONS



EASY LEG TURNER



INFORMATION SHEET

Simple but brilliant, the **Easy Leg Turner** allows a carer to safely raise a person's legs and swivel them into or out of bed – all in a matter of seconds.



EASY LEG TURNER

The **Easy Leg Turner** is designed to assist all single carers to help someone into bed and out of bed. Perfect for supporting the lower limbs whilst raising legs and swivelling in bed or on the sofa. Comfortable to sit on, the person being helped just requires upper sitting balance for a short period of time.

The **Easy Leg Turner** has a Safe Working Load (SWL) of 100kg (15.74 stone) and can be used with a hoist.

Add **Extension Straps**, which carers can attach to the leg lifter to prevent having to bend down when manoeuvring patients/clients. This helps to protect the carers knees and back.



BENEFITS AND FEATURES

BENEFITS

- Rotates really easily
- Excellent support for feet and upper limbs as well
- Easy to remove
- Promotes patient comfort and safety whilst legs are being raised into bed

FEATURES

- Promotes safety, dignity and comfort when getting into/out of bed
- Washable
- One size: 55cm x 200cm
- Suitable for single carers
- Easy to transport, packed in its own carry bag



VISIT THE WEBSITE

QUOTES IN 3 MINUTES

PRODUCT VIDEOS

SECTOR SPECIFIC PRICING

FOR A **5%** DISCOUNT VISIT OUR WEBSITE AND USE THE CODE **5-PERCENT-OFF**

WHY IT'S THE BEST

The **Easy Leg Turner** is an all-in-one solution to getting legs in and out of beds. Ergonomically designed with both the carer and client in mind, the **Easy Leg Turner** is a solution for carers across a wide range of ages and abilities.

The product is both comfy for the client, and designed for the ease of use of the care giver... if required. Additionally, the **Easy Leg Turner** can also be used solo; providing the user has enough upper body strength to support themselves with our extension straps.

Combined with excellent comfort, anti-slip security, and slide-ability on the swivel section, you'll find the **Easy leg Turner** an excellent product to add to your arsenal of products used for getting legs into bed with dignity in mind.



Easy Leg Turner	Size	Product Code
Easy Leg Turner	55cm x 200cm	ELT100200

Wash at 40°C – 74°C
Tumble dry on low heat
Do not iron
Do not use bleach
Do not use fabric conditioner



COMPLIMENTARY PRODUCTS

WEBBING BED LADDERS / ROPE LADDERS


Ideal for maintaining independence and offers support during rehabilitation, enabling patients to sit up and lie down in bed without assistance from a carer.



WRL/10

EXTENSION STRAP – WIPECLEAN

Wipeclean Extension Straps reduce the risk of caregivers over-stretching and over-reaching during supine patient transfers.



LS0575W

BED LOOPS

For patients who struggle to turn in bed, these Bed Loops can be used to securely reposition oneself and safely move to the edge of the bed.



BL-02

TO BUY THE **EASY LEG TURNER**, OR VIEW OUR FULL RANGE, VISIT OUR **WEBSITE** AT

WWW.HOSPITALDIRECT.CO.UK

WHY YOU SHOULD BUY FROM US



Specialising in Patient Moving and Handling across all healthcare sectors, including Hospitals, Community, and Home Care settings; Hospital Direct has been at the forefront of Patient Handling for over 25 years and is the UK's leading manufacturer and supplier of patient handling equipment.

Hospital Direct are internationally recognised as Experts In Patient Handling Solutions, and has built its reputation on product quality, fair pricing, customer service, and fast delivery times. We are also known for having the widest self-manufactured moving and handling product range in the world, and offer professionals access to the 'HD Knowledge Zone' where they can receive training, knowledge, and advice, for situations which require dignity and care. We take our product solutions very seriously, which have all been developed on the back of real issues that healthcare professionals and individuals have required resolutions for.

From helping to get legs into bed in a homecare/community setting, to proning patients in theatres and surgery, Hospital Direct offer a range of standard, exclusive, and bespoke handling solutions which ensures that you and your colleagues can find the right products to fit your needs.

UK: Medical Devices Regulations 2002 (SI 2002 No 618, as amended) (UK MDR 2002).

EU: Class I non-invasive devices registered under or complying with the Medical Devices Directive (MDD) 93/42/EEC (as amended in 2007).

Such a pleasant experience in dealing with Hospital Direct. A company who actually cares.



HOSPITAL



COMMUNITY



RETAIL STORE



EXCLUSIVE PRODUCTS



BUSINESS TRADE AND EXPORT



TRAINING AND DEMOS



Visit our website at
www.hospitaldirect.co.uk

Call: +44 (0)1588 640144

For sales orders: sales@hdigroup.co.uk

For information: info@hdigroup.co.uk

Hospital Direct (Marketing) Ltd
Units 3 & 4, The Green
Clun, Shropshire SY7 8LG



January 2022