




EXPERTS IN PATIENT HANDLING SOLUTIONS



# EASY GUIDE – SATIN SYSTEMS

Our In-bed Systems, such as our SIBA Satin or Flexicare Ranges, help people moving in bed. These patient positioning devices also help carers by enabling them to support loved ones or clients more easily.

## WHAT SATIN IS RIGHT FOR ME?



### 2-WAY BASE SHEET

- Suitable for independent use.
- Use with or without a carer.
- For easier turning in bed (left and right).
- Restricted movement up/down bed when profiled.



### 4-WAY TOP SHEET


- Suitable for use with a care giver.
- Must be used with a base sheet.
- For 4-Directional movement/repositioning.
- Cotton edges act as break on satin when bed is profiled.

## WHAT BASE SHEET DO I NEED? THREE DESIGNS AVAILABLE


### Fitted




### Elastic Straps




### Tuck In



## SATIN PANELS EXPLAINED




MIDI LENGTH SATIN PANEL




FULL LENGTH SATIN PANEL

## EXAMPLES


MIDI (CENTRAL) ON A DOUBLE MATTRESS




MIDI (OFFSET) ON A DOUBLE MATTRESS



FULL LENGTH (OFFSET) ON A DOUBLE MATTRESS



FULL LENGTH (CENTRAL) ON A DOUBLE MATTRESS



## WHAT DESIGN OF TOP SHEET DO I NEED?

1. Top Sheet without Handles
2. Top Sheet with Handles
3. Top Sheet with Locking Loops

## FEATURES AND BENEFITS

The satin we use in all our SIBA Top Sheets and Base Sheets has been specifically woven to be **breathable**, enhancing your patient care and comfort.

Safe Working Load (SWL) – **400kg**.

**Bespoke variants can be made to order.**



# SIBA SATIN SHEET SYSTEM

## BASE SHEETS

**BASE SHEET – SINGLE BED – 90cm WIDE MATTRESS** (Elastic Strap Corner Fitting available – just add ES to the end of the codes)

STMF90200	Fully Fitted Base Sheet – Static and Dynamic Mattress up to 5" depth
STMF90220	Fully Fitted Base Sheet – Static and Dynamic Mattress up to 5" depth cell and extended with SWAB fitted
STMF105200	Fully Fitted Base Sheet – Deep Cell 8" Dynamic Mattress

**BASE SHEET – 4' BED – 120cm WIDE MATTRESS** (Elastic Strap Corner Fitting available – just add ES to the end of the codes)

STMF120200	Fully Fitted Base Sheet – Static and Dynamic Mattress up to 5" depth – Central Satin
STMF120220	Fully Fitted Base Sheet – Static and Dynamic Mattress up to 5" depth cell and extended with SWAB fitted – Central Satin

**BASE SHEET – DOUBLE BED – 135cm WIDE MATTRESS** (Base Sheets with off-set Satin left side or right side) available by request

STMF135200/C	Fully Fitted Base Sheet – Central Satin – 135cm (w) x 200cm (l)
STMF135200/O	Fully Fitted Base Sheet – Offset Satin – 135cm (w) x 200cm (l)

## BASE SHEETS – KING SIZE AND SUPER KING BEDS

STMF150200	Fully Fitted Base Sheet – King Size – full satin area – 150cm (w) x 200cm (l)
STMF180200	Fully Fitted Base Sheet – Super King Size – full satin area – 180cm (w) x 200cm (l)
STMF150200/O	Fully Fitted Base Sheet – King Size – offset satin area – 150cm (w) x 200cm (l)
STMF180200/O	Fully Fitted Base Sheet – Super King Size – offset satin area – 180cm (w) x 200cm (l)

## TOP SHEETS

### 2 WAY TOP DRAW SHEETS – SINGLE BED

ST2W200140	2 Way Top Draw Sheet – Single Bed – 200cm (w) x 140cm (l)
ST2W200200	2 Way Top Draw Sheet – Single Bed – 200cm (w) x 200cm (l)

### 4 WAY TOP DRAW SHEETS WITH OR WITHOUT HANDLES – SINGLE BEDS

ST4W200140H	4 Way Top Draw Sheet with Handles – ¾ Length – 200cm (w) x 140cm (l)
ST4W200140	4 Way Top Draw Sheet without Handles – ¾ Length – 200cm (w) x 140cm (l)
ST4W200200H	4 Way Top Draw Sheet with Handles – Full Length – 200cm (w) x 200cm (l)
ST4W200200	4 Way Top Draw Sheet without Handles – Full Length – 200cm (w) x 200 (l)

### 4 WAY TOP DRAW SHEETS WITH HANDLES – 4' BED

ST4W225140H*	4 Way Top Draw Sheet with Handles – ¾ Length – 225cm (w) x 140cm (l) (single or 4' bed)
ST4W225200H*	4 Way Top Draw Sheet with Handles – Full Length – 225cm (w) x 140cm (l) (single or 4' bed)
ST4W250140H	4 Way Top Draw Sheet with Handles – ¾ Length – 250cm (w) x 140cm (l)
ST4W250200H	4 Way Top Draw Sheet with Handles – Full Length – 250cm (w) x 200cm (l)

\* 4-way sheets for clients that fill the width of a single bed or in a 4' wide bed and have approx 12" space on each side of their body

### 4 WAY TOP DRAW SHEETS WITH OR WITHOUT HANDLES – DOUBLE AND KING-SIZE BED

ST4W250140H	4 Way Top Draw Sheet with Handles – ¾ Length – 250cm (w) x 140cm (l) (4' 6" wide mattress)
ST4W270140	4 Way Top Draw Sheet without Handles – ¾ Length – 270cm (w) x 140cm (l) (4' 6" or 5' wide mattress)
ST4W250200H	4 Way Top Draw Sheet with Handles – Full Length – 250cm (w) x 200cm (l) (4' 6" wide mattress)
ST4W270200	4 Way Top Draw Sheet without Handles – Full Length – 270cm (w) x 200cm (l) (4' 6" or 5' wide mattress)

### 4 WAY TOP DRAW SHEET WITH HANDLES AND LOCKING LOOPS

ST4W200150H/LL/LS2	4 Way Top Draw Sheet with Handles and Locking Loops – ¾ Length – 200cm (w) x 150cm (l)
ST4W200200H/LL/LS2	4 Way Top Draw Sheet with Handles and Locking Loops – Full Length – 200cm (w) x 200cm (l)
ST4W225150H/LL/LS2	4 Way Top Draw Sheet with Handles and Locking Loops – ¾ Length – 225cm (w) x 150cm (l)
ST4W225220H/LL/LS2	4 Way Top Draw Sheet with Handles and Locking Loops – Full Length – 225cm (w) x 220 (l)
ST4W250150H/LL/LS2	4 Way Top Draw Sheet with Handles and Locking Loops – ¾ Length – 250cm (w) x 150cm (l)
ST4W250200H/LL/LS2	4 Way Top Draw Sheet with Handles and Locking Loops – Full Length – 250cm (w) x 200cm (l)

ST4W/LS Locking Straps for above Draw Sheets only

## PILLOWCASE

SPC001

Used for: Prevention and reduction of ear ulcers, helping users move up the bed and moving head while turning. Also use as an aid to prevent patients slipping down in bed.